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
Quantification Methodology for Enzyme Activity Related to Insecticide Resistance in *Aedes aegypti*

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**Quantification Methodology for
Enzyme Activity Related to Insecticide
Resistance in *Aedes aegypti***

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Isabela Reis Montella – IOC/FIOCRUZ - Instituto de Biologia do Exército

Consultor Técnico:

William G Brogdon - Centers for Disease Control and Prevention – CDC/Atlanta

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Fotografias:

Rossana Henriques e Ruy Barbosa Pinto

Impresso no Brasil / Printed in Brazil

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Em inglês: Quantification Methodology for Enzyme Activity Related to Insecticide Resistance in *Aedes Aegypti*

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INDEX

PREFACE	9
ACKNOWLEDGEMENTS	13
1. INTRODUCTION	15
2. ABBREVIATIONS	19
3. EQUIPMENT AND REAGENTS	21
3.1. Equipment	21
3.2. Glassware and plasticware	21
3.3. Miscellaneous	23
3.4. Reagents	23
3.4.1. reagents of general use	23
3.4.2. specific reagents	25
4. PREPARING MOSQUITOES	27
4.1. Freezing	27
4.2. Homogenization	27
4.3. Distribution of the homogenates in the microplates	31
4.3.1. preparing the “mosquito plate”	31
4.3.2. distributing samples in the test microplates	33
5. ENZYMATIC TESTS	35
5.1. Acetylcholinesterase (ACE)	37
5.2. Mixed Function Oxidases (MFO)	39
5.2.1. cytochrome C standard curve	41
5.3. alpha-Esterase (α-EST)	41
5.3.1. alpha-naphthol standard curve	45
5.4. beta-Esterase (β-EST)	45
5.4.1. beta- naphthol standard curve	47
5.5. “PNPA” Esterase (PNPA-EST)	49
5.6. Glutathione-S-transferase (GST)	49
5.7. Total proteins (PTN)	51
5.7.1. protein standard curve – alternative 1	53
5.7.2. protein standard curve – alternative 2	53
5.7.3. distribution of BSA samples for the standard curve	55
6. PREPARING SPECIFIC REAGENTS	57
6.1. Acetylcholinesterase (ACE)	57
6.2. Mixed Function Oxidases (MFO)	59
6.3. alfa-Esterase (α-EST)	61
6.4. beta-Esterase (β-EST)	63

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6.5.	“PNPA” Esterase (PNPA-EST)	65
6.6.	Glutathione-S-transferase (GST)	67
6.7.	Total proteins (PTN)	69
7.	GENERAL BUFFERS AND STOCK SOLUTIONS	71
7.1.	Sodium phosphate buffers	71
7.1.1.	sodium phosphate stocks (mono and bibasic) at 1M	71
7.1.2.	sodium phosphate stock buffers at 1 M	73
7.1.3.	sodium phosphate buffers at 100 and 50 mM	73
7.1.4.	sodium phosphate buffer at 20 mM	73
7.2.	Potassium phosphate buffers	75
7.2.1.	potassium phosphate stocks (mono and bibasic) at 1M	75
7.2.2.	potassium phosphate stock buffers at 1 M	77
7.2.3.	potassium phosphate buffer at 90 mM pH 7.2	77
7.2.4.	potassium phosphate buffer at 100 mM pH 6.5	77
7.3.	Sodium acetate buffer	79
7.3.1.	sodium acetate buffer, stock at 250 mM pH 5.0	79
7.4.	Propoxur (Acetylcholinesterase assay)	79
7.5.	SDS (alpha and beta-Esterases test)	79
8.	COMPARING WHO – CDC – LAFICAVE PROTOCOLS	81
9.	ANALYSIS OF RESULTS	85
9.1.	Acetylcholinesterase (ACE)	85
9.2.	Total proteins (PTN)	85
9.3.	Mixed Function Oxidases (MFO)	87
9.4.	alpha and beta-Esterases (α and β -EST)	89
9.5.	“PNPA” Esterase (PNPA-EST)	89
9.6.	Glutathione-S-transferase (GST)	91
10.	LAFICAVE METHODOLOGY FOR THE PROCESSING AND ANALYSIS OF RESULTS	93
10.1.	File 1 – data transfer from Soft Max to Excel	95
10.2.	File 2 – conversion of absorbance values (except PNPA-EST)	99
10.3.	File 3 – conversion of PNPA-EST absorbance values	109
10.4.	File 4 – total data of a given population	115
11.	INTERPRETATION OF RESULTS	119
11.1.	Criteria for inclusion of evaluations from individual mosquitoes	119
11.2.	Criteria for population classification	121
12.	REFERENCES	127

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
PREFACE


One of the tools employed in the dengue control is the application of insecticides to control its vector. In this context, monitoring the resistance status of *Aedes aegypti* populations, as well as detecting biochemical mechanisms implicated with resistance are very important parameters to define control strategies of this vector.

The publication of this handbook that gives detailed instructions to the quantification, in *Aedes aegypti*, of the activity of enzymes related to insecticide resistance, through methodologies adapted to our conditions, fulfills an important gap and will enable the Laboratories of the Brazilian *Aedes aegypti* Resistance Monitoring Network to accomplish these biochemical assays coordinately and in a better way.

On the other hand, monitoring activities will provide the Health Surveillance Secretary/Brazilian Dengue Control Program General Coordination relevant information to help the rational control of the dengue vector in the country.

With this publication, the Health Surveillance Secretary, together with the Laboratory of Physiology and Control of Arthropod Vectors/Oswaldo Cruz Foundation/Rio de Janeiro, presently National Reference Laboratory of the Brazilian *Aedes aegypti* Insecticide Resistance Monitoring Network, accomplishes its role in the scope of the National Health Surveillance System, contributing to a higher effectiveness in the control actions of an important disease in the context of the epidemiological profile of the country.

Fabiano Geraldo Pimenta Júnior
Health Surveillance Secretary

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
PREFACE

The initial detection and characterization of an instance of insecticide resistance is generally accomplished using some form of bioassay. Bioassays are conducted to answer the question: Will this insecticide control this vector at this place at this time? This is the fundamental practical question that operational control personnel must have answered, preferably before the insecticide is purchased.

Where resistance is found, the question then becomes: What do I do now? It is to answer this question that biochemical and molecular methods were developed for detecting and assessing in individual insects insecticide resistance mechanisms and the potential for cross resistance and its specificity for particular types of insecticides.


Microplate based biochemical resistance assays were initially developed at the Centers for Disease Control in Atlanta in the late 70s. As they were more widely used, a number of modifications in protocol occurred. These were mainly due to the varying filter sets available in commercial micropolates-reading spectrophotometers and to the convenience in particular labs of using different proportions of reagents, more convenient buffers and dilution factors of insect homogenates. In general, two protocols have perhaps been most widely used: the CDC and World Health Organization protocols.

However, the strength of any protocol lies in its adaptability to fulfill the needs of specific control programs. For that reason, modified protocols and original protocols are to be encouraged where they fill an operational need in a country or program. The only caveat is that the assay protocol represent accurately and precisely the biochemistry of insecticide resistance in individual insects.

The LAFICAVE protocol described herein fulfills these requirements, as I have personally verified, and should prove to be a valuable addition to our tool kit of methods for addressing the crucial issue of detecting and describing in practical terms resistance in vectors of disease.

The authors are seeking to provide a practical lab-oriented manual, including what will be a most useful CD rom, that will leave the user with the secure feeling that they know exactly what they are doing when performing these tests. Great care has been exercised in an empirical assessment of every aspect of the methodology. At CDC, we are planning to link this protocol to our resistance training website, further emphasizing that a robust methodology can be tailored to fit the needs of programs.

*William G. Brogdon, Ph.D.
Senior Fellow and Research Entomologist
Entomology Branch
Division of Parasitic Diseases
National Center for Infectious Diseases
Centers for Disease Control and Prevention
Atlanta, Georgia, U.S.A. 30341*

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ACKNOWLEDGEMENTS


Many professionals and friends collaborated to make the publication of the present handbook possible and it would be very difficult without forgetting someone. However, some persons played an important role towards the concretization of this project and must be named. Dr Ima Aparecida Braga, technical assistant of the Brazilian Dengue Control Program/Health Surveillance Department, (PNCD/SVS) encouraged all the steps of both the protocol standardization work and, later, the publication of this volume. Her professional and personal persistence and her friendship were essential to the accomplishment of this project. Equally important has been the support of Dr Giovanini Coelho, coordinator of the Brazilian Dengue Control Program, who, as well as Dr Braga, is a professional extremely engaged with the quality of the public health services in our country, and always takes into account technical and scientific knowledge advances. It is also ought mentioning Dr Fabiano Pimenta, Secretary of the Brazilian Health Surveillance Department, and Dr Jarbas Barbosa, Health Ministry Executive Manager, for their effort to support national abilities. We owe to these four professionals, among many efforts to stimulate the dengue control quality in Brazil, the initiative to bring Dr William G. Brogdon, from the Centers for Disease Control and Prevention, Atlanta, US, which, as an external consultant, personally evaluated all the steps of the protocol and approved the procedures standardized by our team. We would like to thank Dr Brogdon for all his patience and valuable suggestions.

Also relevant to the accomplishment of this work was the Oswaldo Cruz Foundation (Fiocruz), specially the Reference Services and Environment Department, represented by Dr Ary Carvalho de Miranda and his adviser, Ana Beatriz Moraes da Silva, from whom we always received enthusiastic support and an expressive effort in the improvement of the Fiocruz Reference Laboratories and the Oswaldo Cruz Institute, the unit to where we belong, which has been working hard in the interaction between basic and applied research and the qualified public services, aiming at a better public health. Also important is our partnership with the Army Biology Institute, under the command of Colonel Benedito Domicio Ferrari, and previously of General Ivan da Costa Garcez Sobrinho, where we found unrestricted support to the accomplishment of our activities.

Besides Fiocruz and the Brazilian Health Surveillance Department/General Laboratories Coordination (SVS/CGLAB) we received financial support of several programs and research funding agencies, like the National Council for Technological and Scientific Development (CNPq), the Rio de Janeiro Research Funding Agency Carlos Chagas Filho (Faperj), and the Public Health Technological Development Program/Fiocruz (PDTSP).

We would like to mention the English revision of the text and the electronic files made by, respectively, Mitchell Raymond Lishon and Dr Walter Oelemann, the excellent final art work of Rossana Henriques and the drawing of our mascot-mosquito done by Daniel Peixoto. Finally, it should be emphasized that this work would never be finalized without the collaboration of all Arthropod Vectors Physiology and Control Laboratory members, participative in all steps, from the initial decision to define a protocol adapted to our needs until the publication of this handbook.

To these persons, and to all the others that collaborated with us to some extent, our sincere gratitude.

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
1. INTRODUCTION


Chemical insecticides have up now played an important role in the insect disease vectors control. In Public Health, the more widely used insecticides act on the insect central nervous system (CNS) and are classified in four major groups, according to their chemical nature: organochlorines, organophosphates, carbamates and pyrethroids. Insects become resistant to these compounds as a result of either an alteration in their target sites or an increased detoxification efficacy. The latter mechanism is known as metabolic resistance, and the increase in the activity of detoxifying enzymes (Glutathione-S-Transferases, Esterases, Mixed Function Oxidases) contributes to the elimination or deactivation of the insecticides circulating the vector's organism.

In Brazil, as a consequence of the use of insecticides since 1967, dengue vector populations are exposed to an intense selection pressure. In this sense the analysis of both the resistance status of wild populations and the mechanisms involved is of relevance to the control of *Aedes aegypti*. Biochemical methods can be responsible for the early detection of altered specimens, yielding information concerning the resistance mechanisms of a given population as well as a pattern of cross resistance to different insecticides.

This handbook was prepared in order to give detailed instructions to the quantification of the enzyme activity related to insecticide resistance in *A. aegypti*. The work is part of MoReNAa (Brazilian *Aedes aegypti* Insecticide Resistance Monitoring Network) activities, a group of laboratories that act together with the PNCD (Brazilian Dengue Control Program) and give support to the rational control of the dengue vector in the country.

The protocols described in this handbook are the result of modifications in the procedures recommended by the World Health Organization (Hemingway 1998) and by the Centers for Disease Control (Insecticide Resistance Workshop 1998) in order to adapt the methodology to our conditions and needs. These protocols have been utilized with success by our group since 2002.

The main parameters evaluated were: mosquito homogenization volume, nature and concentration of buffers and reactives, size of the aliquots of mosquito homogenate used in each assay and reaction time. We opted to use duplicate samples for all enzymes. We also evaluated the effect of centrifugation on the quantification of Mixed Function Oxidases, since these enzymes are known to be associated with membranes (Montellano 2004). Furthermore, internal controls, corresponding to mosquitoes of the insecticide susceptible reference strain (we use Rockefeller mosquitoes), were included in each assay.

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

We intend to evaluate other parameters, like the differential expression between males and females, throughout development or in adult mosquitoes of varying ages. Currently we work with one-day old adult females. As far as biochemical assays with adults are concerned, CDC methodology utilizes only females while WHO protocol also recommends the use of males. By contrast, WHO protocol works only with young adults, one-day old, while CDC has no special recommendation on this issue. However, temporal differences in the expression of enzymes related to resistance can occur. Differential expression between larvae and adults is also an important parameter to resistance monitoring when it is necessary to compare bioassays and biochemical tests.

We defined protocols to quantify Acetylcholinesterase, the target site in the CNS of organophosphates and carbamates, and the main enzymes related to metabolic resistance. In conjunction, these enzymes can determine resistance to the major insecticide classes employed in Public Health. We describe the quantification of Esterases with three different substrates. All the protocols are based on the quantification of the enzyme activity. The exception is the Mixed Function Oxidases (MFO) assay, an indirect measure of the enzyme content, since it quantifies the heme group, its prosthetic group (see [Section 9.3](#)).


To facilitate handling, the text is organized in separate sections, from freezing of mosquitoes to interpretation of the results. Details of preparation of buffers and reagents can also be found in separate sections. Finally, four Excel files accompany this handbook. These files enable, in a semi-automatic way, to compare replicas of each sample, to calculate the enzymatic activity of each mosquito, to correct values obtained by the total protein amount of each mosquito and to elaborate graphs with the activity profile of each tested population.

We hope this work contributes to the *Aedes aegypti* insecticide resistance monitoring that is presently being performed in the country, stimulating other laboratories to perform biochemical assays. It could also serve as an initial guide to the evaluation of resistance mechanisms in other insect species. Initially planned as a handbook to be used by the MoReNAa network, we were later encouraged to present this material to a wider public; this is the reason for the English version.

2. ABREVIATURAS

ACE	Acetilcolinesterase
AChE	atividade total de Acetilcolinesterase
AChI	atividade de Acetilcolinesterase na presença de inibidor
ACTH	iodeto de acetiltiocolina
BSA	albumina sérica bovina
cit	citocromo
CDNB	1-cloro 2,4 dinitrobenzeno
DIA	orto dianisidina, tetrazotizada
DTNB	ácido 5,5 ditiobis (2-nitro-benzóico)
EST	Esterase
GSH	L-glutationa reduzida
GST	Glutationa-S-transferase
H₂O₂	peróxido de hidrogênio
K fosf	fosfato de potássio
MFO	Oxidases de Função Mista
Na acet	acetato de sódio
PNPA	4-nitrofenil acetato
PTN	proteínas
TMBZ	3,3,5,5 dihidrocloro de tetrametil benzidina

2. ABBREVIATIONS

ACE	Acetylcholinesterase
AChE	total activity of Acetylcholinesterase
AChI	Acetylcholinesterase activity in the presence of an inhibitor
ACTH	acetylthiocholine iodide
BSA	bovine serum albumin
cit	cytochrome
CDNB	1-chloro-2,4-dinitrobenzene
DIA	o-dianisidine, tetrazotized
DTNB	5,5'-Dithiobis (2-nitrobenzoic acid)
EST	Esterase
GSH	glutathione
GST	Glutathione-S-transferase
H₂O₂	hydrogen peroxide
K phosph	potassium phosphate
MFO	Mixed Function Oxidases
Na acet	sodium acetate
PNPA	p-nitrophenyl acetate
PTN	proteins
TMBZ	3,3',5,5'-tetramethyl-benzidine dihydrochloride

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3. EQUIPMENT AND REAGENTS

3.1. Equipment

- Ice machine (flaked ice)
- Water purification system (reagent grade water)
- Microcentrifuge for Eppendorf tubes
- Microplates spectrophotometer (Versa-max, Molecular Devices)
- Single channel pipettes, P20, P200 and P1000, variable volumes
- Multichannel pipette (12 channels), 5 to 50 μ l
- Multichannel pipette (12 channels), 30 to 300 μ l or 50 to 300 μ l
- pH Meter
- Analytical balance
- Refrigerator
- Freezer -20°C
- Freezer -80°C
- Magnetic stirrer

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3.3. Miscellaneous

- Filtration system
- Homogenizer

We use an electric homogenizer, with pedal. A battery homogenizer also can be employed. The same pestle can be used in both cases.

Battery homogenizer.
Facing page:
electric homogenizer,
with pedal

3.4 Reagents¹

3.4.1. reagents of general use

reagent	manufacturer	cat n°	quantity
Sodium acetate, anhydrous	Sigma	S-2889	1 Kg
Monobasic potassium phosphate, anhydrous	Sigma	P-5379	500 g
Bibasic potassium phosphate, anhydrous	Sigma	P-3786	
Bibasic potassium phosphate, tri-hidratad	Sigma	P-5504	
Monobasic sodium phosphate, anhydrous	Sigma	S-5011	1 Kg
	Sigma	S-0751	500 g
Bibasic sodium phosphate, anhydrous	Sigma	S-5136	1 Kg

Phosphate buffers can be prepared with anhydrous or hydrated phosphates, providing the molecular weight of the original reagents is considered.

All the reagents listed above can be stored at room temperature.

¹ Use of trademark names does not imply endorsement by the Health Ministry but is intended only to assist in identification of a specific product. Use of reactivates from other manufacturers should be preceded by evaluation tests.

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(*) light sensitive reagent. Store protected from light.

(**) Although storage of these reagents at room temperature is originally recommended, in the laboratory we keep them at 4°C, because of our tropical climate.

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4. PREPARING MOSQUITOES

4.1. Freezing

- a) Use one-day old (0 - 24 hours after emergence) adult females, NON blood fed².
- b) Females from the populations to be tested are placed still alive inside Eppendorf tubes in pools of 40-50. Females of an insecticide susceptible reference strain, used as internal controls of the assays, are stored in groups of seven (7).
- c) Tubes containing mosquitoes are immediately frozen (*females still alive*) at -70°C ³.
- d) When necessary, remove tubes from the freezer *immediately before use* (Section 4.2, c). Maintain the tubes iced. It is not recommended to re-freeze the remaining mosquitoes.

4.2. Homogenization

- a) Distribute 45 numbered Eppendorf tubes in a basin containing ice.
- b) Add 30 μl of water to each tube. Use preferably reagent grade water (ex: Milli-Q)⁴.
- c) Remove one tube containing 45-50 mosquitoes to be tested and one tube with the insecticide susceptible mosquitoes from the freezer.
- d) Spread mosquitoes over a plastic basin (or a Petri dish), directly placed over ice.
- e) Put one mosquito in each numbered tube, 1-40 corresponding to mosquitoes from the population to be tested and 41-45, susceptible mosquitoes, from the reference strain.

² According to the WHO handbook (Hemingway 1998), it is not recommended to use mosquitoes that were directly exposed to pesticides.

³ According to the WHO handbook (Hemingway 1998), dead insects kept at room temperature for more than some minutes have enzyme levels seriously reduced.

⁴ Do not store neither autoclave Milli-Q water. Reverse osmosis purified water can also be used.

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- f) Homogenize each mosquito^{5,6}, following the numbers defined for the tubes. *Keep tubes iced during the whole homogenization procedure*⁷.
- g) Add 270 µl of water to each tube. DO NOT pipette on the tube wall: if you pipette in the middle of the tube, mixing will be favored. Respect the numbers defined for the tubes (1-45). In this way, when the last tube is attained, the cell debris from the first samples will already have sunk to the bottom of the tube. You can now start the distribution of aliquots to the Acetylcholinesterase (ACE) and Mixed Function Oxidases (MFO) assays (**Sections 5.1 and 5.2**).
- h) Before centrifugation, remove four aliquots, 25 µl each, to quantify ACE; each pair of aliquots in one separate 96-well microplate (see **Section 4.3.2**)⁸.
- i) Before centrifugation, transfer two additional aliquots, 20 µl each, to one 96-well microplate, in order to quantify MFO⁸.
- j) Centrifuge the remaining homogenate (~160 µl) at 12.000 g for 60 seconds^{9,10}.
- k) Remove tubes from the centrifuge and transfer to the basin containing ice. Keep tubes on ice while you distribute the homogenates, in order to reduce proteolysis¹¹.
- l) Transfer to a new 96-well microplate (the “mosquito plate”, **Section 4.3.1**), about 130-150 µl of the supernatant from each homogenate¹². From this point on, a multi-channel pipette can be used to distribute the aliquots.

Facing page:
homogenization
detail of an individual
mosquito, with an electric
homogenizer

⁵ Automatic homogenizers are very efficient. We use an adapted bench equipment, with pedal. Portable homogenizers, battery-operated, are also indicated (**Section 3.3**).

⁶ Each mosquito is homogenized during 5-10 seconds, approximately. After homogenizing each mosquito, the pestle is washed in a 100 ml capacity Beaker, containing 50 ml water. Water in the Beaker is not changed during all the procedure. After each washing, the pestle is dried with a Kleenex.

⁷ *Once you started the homogenization, it is important to proceed quickly.*

⁸ Avoid to resuspend and to remove the decanted material at the bottom of the tube.

⁹ According to the WHO handbook (**Hemingway 1998**), the centrifugation step is important mainly to the GST dosage (particles can interfere with absorption at 340 nm, see **Section 5.6**).

¹⁰ Ideally use a refrigerated centrifuge. If only non-refrigerated centrifuges are available, do not centrifuge for more than 2 minutes. We use a microcentrifuge for Eppendorf tubes, at maximum speed.

¹¹ According to the WHO handbook (**Hemingway 1998**), freezing at -20°C immediately after the homogenization procedure keeps enzyme activity relatively constant for several hours. However, in our experience, it is important to perform the test as quickly as possible.

¹² A total of 110 µl of each homogenate will be needed in the remaining quantifications (besides Acetylcholinesterase and MFO, already discriminated, items h and i).

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4.3. Distribution of the homogenates in the microplates

4.3.1. preparing the “mosquito plate”

To speed up the tests, 130-150 μ l of the supernatant from the homogenates of individual mosquitoes are transferred from the Eppendorf tubes to a 96-well microplate (always on the ice), in the following order¹³:

	1	2	3	4	5	6	7	8	9	10	11	12
A	m1	m2	m3	m4	m5	m6	m7	m8	m9	m10	m11	m12
B	m13	m14	m15	m16	m17	m18	m19	m20	m21	m22	m23	m24
C	m25	m26	m27	m28	m29	m30	m31	m32	m33	m34	m35	m36
D	m37	m38	m39	m40	R41	R42	R43	R44	R45			
E												
F												
G												
H												

m = mosquito from the population to be tested

R = mosquito from the insecticide susceptible reference strain

Facing page: transfer of mosquito homogenate supernatants to the “mosquito plate”

¹³ This distribution enables the use of the same microplate in two tests: lines A-D and lines E-H.

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4.3.2. distributing samples in the test microplates

Prepare five (5) different microplates, containing the following amounts of mosquito homogenate, in duplicate samples:

enzyme	aliquots
GST	2 X 15 µl
alpha-Esterase	2 X 10 µl
beta-Esterase	2 X 10 µl
“PNPA” Esterase	2 X 10 µl
Total proteins	2 X 10 µl

Distribute the homogenate aliquots according to the model below:

	1	2	3	4	5	6	7	8	9	10	11	12
A	m1	m2	m3	m4	m5	m6	m7	m8	m9	m10	m11	m12
B	m1	m2	m3	m4	m5	m6	m7	m8	m9	m10	m11	m12
C	m13	m14	m15	m16	m17	m18	m19	m20	m21	m22	m23	m24
D	m13	m14	m15	m16	m17	m18	m19	m20	m21	m22	m23	m24
E	m25	m26	m27	m28	m29	m30	m31	m32	m33	m34	m35	m36
F	m25	m26	m27	m28	m29	m30	m31	m32	m33	m34	m35	m36
G	m37	m38	m39	m40	R41	R42	R43	R44	R45	b	b	b
H	m37	m38	m39	m40	R41	R42	R43	R44	R45	c	c	c

m = mosquito from the population to be tested

R = mosquito from the insecticide susceptible reference strain

b = blank (negative control)

c = positive control¹⁴

Facing page: before centrifugation, mosquito homogenate aliquots are transferred to 96-well microplates in order to quantify ACE and MFO

¹⁴ When no positive control is available, the 6 wells (G10 to G12 and H10 to H12) are used as blanks (see Section 5).

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5. ENZYMATIC TESTS

In order to optimize the procedure, each test is always performed simultaneously by two persons. We DO NOT recommend that a complete test (45 mosquitoes, all enzymes simultaneously) be done by just one person.

The test for each enzyme has a different protocol and a different reaction time. In the laboratory we defined an optimized sequence for all the reactions:

- 1) Mixed Function Oxidases¹⁵
- 2) Acetylcholinesterase
- 3) alpha-Esterase¹⁵
- 4) beta-Esterase¹⁵
- 5) total proteins
- 6) protein standard curve¹⁶
- 7) “PNPA” Esterase
- 8) Glutathione-S-transferase

This sequence of reactions is not obligatory, excluding MFO and ACE that are the first enzymes to be quantified. Each laboratory should define the sequence of reactions most convenient for its routine.

¹⁵ In addition to the quantification of this enzyme, it is necessary to do standard curves with the substrates of MFO and alpha- and beta-Esterases. These curves are necessary to the conversion of absorbance values in the specific activity of each enzyme. In our routine, we opted to construct new standard curves each time a new batch of tests with field populations starts (10-20 populations). *We take into account the mean of at least three standard curves.*

¹⁶ A new protein standard curve is done at each test. It can be done at any moment during the test. However, it is recommended that it is performed together with the protein dosage of the mosquito homogenates.

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
Resulting color of ACE assay, in the AChE (left panel) and AChI (right panel) plates, respectively with and without propoxur. B: blank. Facing page: mosquito homogenate distribution in the microplates for ACE test

5.1. Acetylcholinesterase (ACE)

The table below discriminates the solutions necessary for this test. More details related to the preparation of solutions are shown in [Section 6.1](#).

n°	solution	quantity per test		observations
		quantity used	quantity prepared	
1	Triton/Na phosphate	29 ml	500 ml	in stock
2	DTNB/Na phosphate	2 ml	3 ml	freshly prepared
3	AChE	2.5 ml	3 ml	
4	AChI	2.5 ml	3 ml	

- Distribute 25 μ l of the homogenates, *before centrifugation*, in duplicate, in two microplates. One plate is named “AChE” and the other “AChI”. In the AChE plate the total activity of Acetylcholinesterase is quantified. In plate AChI the remaining activity of this enzyme, in the presence of an inhibitor, is measured. Propoxur is the inhibitor utilized, carbamate indicated by the WHO handbook ([Hemingway 1998](#)).
- In the “blank” wells (6 wells: G10 to G12 and H10 to H12), add 25 μ l of water.
- Add to each well 145 μ l of Triton / Na phosphate.
- Add to each well 10 μ l DTNB / Na phosphate.
- Add to each well of the plate AChE 25 μ l of the solution AChE.
- Add to each well of the plate AChI 25 μ l of the solution AChI.
(the final volume in each well is 205 μ l)
- Read the absorbance of both plates (AChE and AChI), at 405 nm after 1 hour (“end point”).

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

Resulting color of MFO assay. B: blank; C: positive control. *Facing page*: mosquito homogenate distribution for MFO test (left panel) and distribution of reagents using a multichannel pipette (right panel)

5.2. Mixed Function Oxidases (MFO)

This test does not measure the enzymatic activity. It quantifies the heme group that, in non blood-fed insects, is mainly associated with cytochrome P450 (see [Section 9.3](#) for details).

The table below discriminates the solutions necessary for this test. More details related to the preparation of solutions are shown in [Section 6.2](#).

n°	solution	quantity per test		observations
		quantity used	quantity prepared	
1	90 mM K phosph, pH 7.2	6 ml	500 ml	in stock
2	0.01 mg/ml cytochrome C	60 µl	5 ml	
3	TMBZ/Na acet	20 ml	24 ml	freshly prepared
4	3% H ₂ O ₂	2.5 ml	3 ml	

- Distribute 20 µl of the homogenates, before centrifugation, in duplicate, in one microplate.
 - In the “blank” wells (3 wells: G10 to G12), add 20 µl of K phosph.
 - In the “positive control” wells (3 wells: H10 to H12), add 20 µl of the solution cytochrome C (total per test: 0.2 µg cytochrome C)¹⁷.
 - Add to each well 60 µl K phosph buffer.
 - Add to each well 200 µl of the working solution: TMBZ/Na acet.
 - Add to each well 25 µl of 3% H₂O₂ (*).
- (the final volume in each well is 305 µl)
- Incubate at room temperature and protected from light during 90 minutes.
 - Read the absorbance at 650 nm (“end point”).

¹⁷ The concentration of cytochrome C used in the positive control was calibrated according to the standard curve and it is included in its linearity range. It is recommended to frequently observe the absorbance value of the positive control, generally very stable among successive tests. It is recommended to repeat the standard curve if the values observed for the cytochrome C are very different from the usual ones.

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5.2.1. cytochrome C standard curve

In our laboratory routine, we opted to devise a new standard curve for every new battery of tests with field populations. We consider the mean of at least three standard curves as a basis to calculate the angular coefficient that is utilized in the conversion of absorbance to μg substrate (Section 9.3).

Add, in duplicate, in the same line of a microplate, the following quantities of 0.01 mg/ml cytochrome C:

μg cit C	μl cit C (0.01 mg/ml)	μl K phosph
0.0	–	80
0.01	1	79
0.02	2	78
0.05	5	75
0.1	10	70
0.2	20	60

Proceed exactly as to the quantification of mosquito samples – Section 5.1, e – h.

5.3. alpha-Esterase (α -EST)

The Table below discriminates the solutions necessary for this test. More details related to the preparation of solutions are shown in Section 6.3.

n°	solution	quantity per test		observations
		quantity used	quantity prepared	
1	alpha-naphthol (positive control)	30 μl	5 ml	in stock
2	alpha-naphthyl acetate/Na phosphate	20 ml	25 ml	freshly prepared
3	Fast Blue	10 ml	15 ml	

- Distribute 10 μl of the homogenates, *after centrifugation*, in duplicate, in one microplate.
- In the “blank” wells (3 wells: G10 to G12) add 10 μl of water.

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
- c) In the “positive control” wells (3 wells: H10 to H12), add 10 μl of the solution alpha-naphthol at 0.5 $\mu\text{g}/\mu\text{l}$ (~ 3.5 nmoles/ μl) in each well (total: 5 μg or ~ 35 nmoles)^{18,19}.
- d) Add to each well 200 μl of alpha-naphthyl acetate/Na phosphate.
- e) Incubate at room temperature during 15 minutes.
- f) Add to each well 50 μl of Fast Blue. Prepare this solution in quantity enough to the dosage of alpha and beta Esterases (5 ml are used in each test). The solution contains a detergent, and foam is produced during its preparation. For this reason it is recommended to prepare a volume higher than the volume used. Agitate by inversion (do not use the vortex). Pipette with care.
(the final volume in each well is 260 μl)
- g) Incubate at room temperature during 5 minutes.
- h) Read the absorbance at 570 nm (end point).

Resulting color of α -EST assay. B: blank; C: positive control. Facing page: mosquito homogenate distribution for α -EST and β -EST tests

¹⁸ The positive controls, alpha- and beta-naphthol, are stable. However, it is convenient to freeze 50 μl aliquots (-20°C) and thaw it just once.

¹⁹ The concentration of alpha- and beta-naphthol used in the positive control (0.5 $\mu\text{g}/\mu\text{l}$) was calibrated according to the standard curve and it is included in its linearity range. It is recommended to frequently observe the absorbance value of the positive control, generally very stable among successive tests. It is recommended to repeat the standard curve if the values observed for the positive controls are very different from the usual ones.

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
5.3.1. alpha-naphthol standard curve

Add, in duplicate, in the same line of a microplate, the following quantities of 0.5 mg/ml alpha-naphthol:

μg alpha-naphthol	μl alpha-naphthol (0.5 mg/ml)	μl H ₂ O
0	–	10
1	2	8
2	4	6
3	6	4
4	8	2
5	10	–

Proceed exactly as to the quantification of mosquito samples – [Section 5.3, d – h](#).

5.4. beta-Esterase (β -EST)

The table below discriminates the solutions necessary for this test. More details related to the preparation of solutions are shown in [Section 6.4](#).

n°	solution	quantity per test		observations
		quantity used	quantity prepared	
1	beta-naphthol (positive control)	30 μl	5 ml	in stock
2	beta-naphthyl acetate/Na phosphate	20 ml	25 ml	freshly prepared
3	Fast Blue	10 ml	15 ml	

- Distribute 10 μl of the homogenates, *after centrifugation*, in duplicate, in one microplate.
- In the “blank” wells (3 wells: G10 to G12) add 10 μl of water.

- g) Esperar 5 minutos à temperatura ambiente.
 h) Fazer a leitura (*end point*), a 570 nm.

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- c) In the “positive control” wells (3 wells: H10 to H12), add 10 μl of the solution beta-naphthol at 0.5 $\mu\text{g}/\mu\text{l}$ (~ 3.5 nmoles/ μl) in each well (total: 5 μg or ~ 35 nmoles).
- d) Add to each well 200 μl of beta naphthyl acetate/Na phosphate.
- e) Incubate at room temperature during 15 minutes.
- f) Add to each well 50 μl of Fast Blue. As mentioned in the previous section, this solution is prepared in quantity enough to the dosage of alpha and beta-Esterases. Pipette with care.
(the final volume in each well is 260 μl)
- g) Incubate at room temperature during 5 minutes.
- h) Read the absorbance at 570 nm (end point).

5.4.1. beta-naphthol standard curve

Add, in duplicate, in the same line of a microplate, the following quantities of 0.5 mg/ml beta-naphthol:

μg beta-naphthol	μl beta-naphthol (0,5 mg/ml)	μl H ₂ O
0	–	10
1	2	8
2	4	6
3	6	4
4	8	2
5	10	–

Proceed exactly as to the quantification of mosquito samples – [Section 5.4, d – h](#).

Facing page: Addition of specific Esterase reagents in 96-well test microplates

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

ROCK

The GST assay is read with a UV filter (340nm) since its reaction product is not colored. Facing page: Resulting color of PNPA-EST assay. B: blank

5.5. "PNPA" Esterase (PNPA-EST)

The table below discriminates the solutions necessary for this test. More details related to the preparation of solutions are shown in [Section 6.5](#).

n°	solution	quantity per test		observations
		quantity used	quantity prepared	
1	PNPA/Na phosphate	20 ml	25 ml	freshly prepared

- Distribute 10 μ l of the homogenates, *after centrifugation*, in duplicate, in one microplate.
- In the "blank" wells (6 wells: G10 to G12 and H10 to H12), add 10 μ l of water.
- Add to each well 200 μ l of the working solution: PNPA/Na phosphate.
(the final volume in each well is 210 μ l)
- Read the absorbance at 405 nm, at 15-second intervals, during 2 minutes.
According to the WHO handbook ([Hemingway 1998](#)), *this assay must not be read as "end point"*.

5.6. Glutathione-S-transferase (GST)

The table below discriminates the solutions necessary for this test. More details related to the preparation of solutions are shown in [Section 6.6](#).

n°	solution	quantity per test		observations
		quantity used	quantity prepared	
1	GSH/CDNB	20 ml	21 ml	freshly prepared

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
Resulting color of PTN assay. B: blank; C: positive control. *Facing page:* Addition of protein reactive to quantify total proteins in mosquito homogenates

- Distribute 15 μl of the homogenates, *after centrifugation*, in duplicate, in one microplate.
- In the “blank” wells (6 wells: G10 to G12 and H10 to H12) add 15 μl of water.
- Add to each well 195 μl of the working solution: GSH/CDNB.
(*the final volume in each well is 210 μl*)
- Read the absorbance at 340 nm, at one-minute intervals, during 20 minutes.

5.7. Total proteins (PTN)

The table below discriminates the solutions necessary for this test. More details related to the preparation of solutions are shown in [Section 6.7](#).

n°	solution	quantity per test		observations
		quantity used	quantity prepared	
1	1 $\mu\text{g}/\mu\text{l}$ BSA	30 μl	1 ml	in stock
2	Bio-Rad reactive	30 ml	40 ml (**)	freshly prepared

(**) part of this quantity is used in the BSA standard curve: 3.6 ml ([Section 5.7.1](#), alternative 1) or 7.2 ml ([Section 5.7.2](#), alternative 2).

- Distribute 10 μl of the homogenates, *after centrifugation*, in duplicate, in one microplate.
- In the “blank” wells (3 wells: G10 to G12) add 10 μl of water.
- In the “positive control” wells (3 wells: H10 to H12), add 10 μl of 1 $\mu\text{g}/\mu\text{l}$ BSA (total:10 μg).
- Add to each well 300 μl of the Bio-Rad reactive at 1:5.
(*the final volume in each well is 310 μl*)
- Read the absorbance 3 -5 minutes later, at 620 nm.

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5.7.1. protein standard curve – alternative 1

This standard curve is routinely used. Add, in duplicate, in the same row of a microplate, as indicated in [Section 5.7.3](#), the following quantities of BSA:

$\mu\text{g BSA}$	$\mu\text{l BSA (*)}$	$\mu\text{l H}_2\text{O}$
0	0	10
2.5	2.5	7.5
5	5.0	5
7.5	7.5	2.5
10	10.0	0
15	7.5	2.5

(*) A 1 $\mu\text{g}/\mu\text{l}$ BSA stock solution is used in all cases, except in the point "15 μg ", when a 2 $\mu\text{g}/\mu\text{l}$ BSA solution is used.

Proceed as to the quantification of mosquito samples – [Section 5.7, d, e](#).

5.7.2. protein standard curve – alternative 2

This standard curve can be used when a higher level of detail is needed. Add, in duplicate, in two rows of a microplate, as indicated in [Section 5.7.3](#), the following quantities of BSA:

$\mu\text{g BSA}$	$\mu\text{l BSA (*)}$	$\mu\text{l H}_2\text{O}$
0	0.0	10
1	1.0	9
2	2.0	8
3	3.0	7
4	4.0	6
5	5.0	5
6	6.0	4
7	7.0	3
8	8.0	2
9	9.0	1
10	10.0	0
15	7.5	2.5

(*) A 1 $\mu\text{g}/\mu\text{l}$ BSA stock solution is used in all cases, except in the point "15 μg ", when a 2 $\mu\text{g}/\mu\text{l}$ BSA solution is used.

Facing page: Example of absorbance reading in the spectrophotometer

Proceed as to the quantification of mosquito samples - [Section 5.7, d, e](#).

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6 PREPARING SPECIFIC REAGENTS

6.1. Acetylcholinesterase (ACE)

n°	solution	quantity per test		observations
		quantity used	quantity prepared	
1	Triton/Na phosph	29 ml	500 ml	in stock
2	DTNB/Na phosph	2 ml	3 ml	freshly prepared
3	AChE	2.5 ml	3 ml	
4	AChI	2.5 ml	3 ml	

1) **Triton/Na phosph** – 1% Triton X-100 in 100 mM sodium phosphate buffer pH 7.8

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Stock, 500 ml

- 5 ml 100% Triton X-100
- 50 ml 1 M sodium phosphate buffer pH 7.8 (details of preparation in [Section 7.1.2](#))
- water to 500 ml

2) **DTNB/Na phosph** – 10 mM DTNB in 100 mM sodium phosphate buffer pH 7.0

This solution must be prepared *just before use* (4 hours maximum).

Working solution, 3 ml

- 0.012 g (0.01189 g) DTNB (MW 396.3)²⁰
- 3 ml 100 mM sodium phosphate buffer pH 7.0 (details of preparation in [Section 7.1.3](#))

3, 4) **AChE/AChI** – 10 mM acetylthiocholine iodide, in water, in the absence (*AChE*) or in the presence (*AChI*) of propoxur.

Both solutions must be freshly prepared (*maximum 4 hours before use*).

In order to have 3 ml of each solution:

Working solution, 6 ml

- 0.0174 g acetylthiocholine iodide (MW 289.2)²¹
- 6 ml water
- reserve 3 ml: *AChE*
- to the remaining 3 ml, add 6 µl 0.1 M propoxur in acetone: *AChI* (details of propoxur preparation in [Section 7.4](#))

²⁰ We keep several 0.012 g DTNB aliquots in Eppendorf tubes, at 4°C.

²¹ We keep several 0.0174 g acetylthiocholine aliquots in Eppendorf tubes, at -20°C.

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a) 0,2% TMBZ em metanol (6 ml por ensaio)

6.2. Mixed Function Oxidases (MFO)

n°	solution	quantity per test		observations
		quantity used	quantity prepared	
1	K phosph	6 ml	500 ml	in stock
2	cytochrome C	60 µl	5 ml	
3	TMBZ/Na acet	20 ml	24 ml	freshly prepared
4	3% H ₂ O ₂	2.5 ml	3 ml	

1) **potassium phosphate buffer (K phosph)** - 90 mM potassium phosphate buffer
pH 7.2

This solution can be stored for several months, closed at 4°C.

Stock, 500 ml

- details of preparation in [Section 7.2.3](#)

2) **cytochrome C** – 0.01 mg/ml in 250 mM sodium acetate pH 5.0

This solution is stored at -20°C, in 80-100 µl aliquots. *Each aliquot is used only once, and the remainder is discarded.*

Stock, 5 ml

- weight 0.5 mg cytochrome C

- add 5 ml 250 mM sodium acetate buffer pH 5.0 ([Section 7.3](#)). Final concentration:
0.1 mg/ml

- dilute this solution to a tenth: 0.5 ml stock + 4.5 ml 250 mM sodium acetate pH 5.0

- this amount is enough for 50 aliquots of 100 µl each

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Working solution

- 0.012 g TMBZ²²
- 6 ml methanol

This solution must be freshly prepared.

b) 250 mM sodium acetate buffer pH 5.0

Stock

- details of preparation in [Section 7.3.1](#)

4) **hydrogen peroxide** – 3% H₂O₂

This solution must be freshly prepared.

Working solution

- 300 µl 30% H₂O₂
- 2.7 ml water

6.3. alfa-Esterase (α -EST)

n°	solution	quantity per test		observations
		quantity used	quantity prepared	
1	alpha-naphthol (positive control)	30 µl	5 ml	in stock
2	alpha-naphthyl acet/Na phosph	20 ml	25 ml	freshly prepared
3	Fast Blue	5 ml	15 ml	

1) **alpha-naphthol** – 0.5 mg/ml alpha-naphthol (~ 3.5 mmoles/ml) (MW: 144.2)

This solution is stored at -20°C, in 100 µl aliquots. Aliquots may be used more than once.

Stock, 5 ml

- weigh 2.5 mg de alpha-naphthol
- add 5 ml acetone
- this amount is enough for 50 aliquots of 100 µl each

2) **alpha-naphthyl acetate/sodium phosphate (alpha-naphthyl acet/Na phosph)** – 0.3 mM alpha-naphthyl acetate in 20 mM sodium phosphate buffer pH 7.2

This solution must be prepared *not more than 1-2 hours before use*.

²² We keep several 0.012 g TMBZ aliquots in Eppendorf tubes, at 4°C.

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Working solution

- 250 μ l 30 mM alpha-naphthyl-acetate (*a*)
- 24.75 ml 20 mM sodium phosphate buffer pH 7.2 (*b*)

This solution is prepared in two steps:

- a*) 30 mM alpha-naphthyl acetate in acetone

Stock

- 0.028 g alpha-naphthyl acetate (MW 186.2)
- 5 ml acetone

This solution can be stored for several months, closed, at 4°C.

- b*) 20 mM sodium phosphate buffer pH 7.2

- details of preparation in [Section 7.1.4](#)

- 3) **Fast Blue** – 0.3% Fast Blue B in 3.5% SDS

This solution must be prepared not more than 1-2 hours before use.

Working solution

- weigh 45 mg Fast Blue B (0.045 g)²³
- add 4.5 ml water
- *After solubilization*, add 10.5 ml 5% SDS ([Section 7.5](#))

6.4. beta-Esterase (β -EST)

n°	solution	quantity per test		observations
		quantity used	quantity prepared	
1	beta- naphthol (positive control)	30 μ l	5 ml	in stock
2	beta- naphthyl acet/Na phosph	20 ml	25 ml	freshly prepared
3	Fast Blue	5 ml	15 ml	

- 1) **beta-naphthol** - 0.5 mg/ml beta-naphthol (~ 3.5 mmoles/ml)

This solution is stored at -20°C, in 100 μ l aliquots. Aliquots can be used more than once.

Stock, 5 ml

- weigh 2.5 mg de beta-naphthol
- add 5 ml acetone
- this amount is enough for 50 aliquots of 100 μ l each

²³ We keep several 45 mg (0.045 g) Fast Blue aliquots in Eppendorf tubes, at 4°C.

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2) **beta-naphthyl acetate/sodium phosphate (beta- naphthyl acet/Na phosph)** – 0.3 mM beta-naphthyl acetate in 20 mM sodium phosphate buffer pH 7.2

This solution must be prepared not more than 1-2 hours before use.

Working solution

- 250 µl 30 mM beta-naphthyl acetate (a)
- 24.75 ml 20 mM sodium phosphate buffer pH 7.2 (b)

This solution is prepared in two steps:

- a) 30 mM beta-naphthyl acetate in acetone

Stock

- 0.028 g beta-naphthyl acetate (MW 186.2)
- 5 ml acetone

This solution can be stored for several months, closed, at 4°C.

- b) 20 mM sodium phosphate buffer pH 7.2

- details of preparation in [Section 7.1.4](#)

- 3) **Fast Blue** – Preparation of this solution, in quantity enough to the alpha and beta-Esterases assays, is indicated in the previous section.

6.5. "PNPA" Esterase (PNPA-EST)

n°	solution	quantity per test		observations
		quantity used	quantity prepared	
1	PNPA/Na phosph	20 ml	25 ml	freshly prepared

- 1) **p-nitrophenyl acetate/sodium phosphate (PNPA/Na phosph)** – 1 mM PNPA in 50 mM sodium phosphate buffer pH 7.4

This solution must be prepared immediately before use (in general our preparation *does not exceeded 5' prior to usage*).

If the solution presents an intense yellow color, discard.

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This solution is prepared in two steps:

a) 100 mM PNPA in acetonitrile

Working solution

- 0.01815 g PNPA (MW 181.15)²⁴

- 1 ml acetonitrile²⁵

According to the WHO handbook ([Hemingway 1998](#)), this solution can be stored during 1-2 months, closed, at 4°C. *In the laboratory we discard after one week.*

Acetonitrile is highly toxic. Attention!!

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1) **Reduced glutathione/chloro-dinitrobenzene (GSH/CDNB)** – 9.5 mM GSH in 1 mM CDBN

The WHO handbook ([Hemingway 1998](#)) recommends preparing this solution 1-2 hours before use. *We prepare at most 10-15 minutes before.*

Working solution

- 20 ml 10 mM GSH (*a*)

- 1 ml 21 mM CDBN (*b*)

This solution is prepared in two steps:

a) 10 mM GSH (MW 307.3) in potassium phosphate buffer

²⁴ This solution is enough to four (4) tests. The solution described in (a) is kept at 4°C. In addition, several 0.01815 g PNPA aliquots can be stored in Eppendorf tubes at -20°C.

²⁵ We already had problems derived from the quality of acetonitrile. We presently we use Sigma's product that works well.

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Working solution

- 0.0615 g GSH²⁶
- 20 ml 100 mM potassium phosphate buffer pH 6.5 (see [Section 7.2.4](#))

b) 21 mM CDNB (MW 202.6) in methanol

Working solution

- 0.0042 g CDNB²⁴
- 1 ml methanol

6.7. Total proteins (PTN)

n°	solution	quantity per test		observations
		quantity used	quantity prepared	
1	BSA (control)	30 µl	1 ml	in stock
2	Bio-Rad reactive	30 ml	40 ml (**)	freshly prepared

1) **Bovine Serum Albumin (BSA)** – 1 mg/ml in water

Vials with 1 ml of the standard BSA solution, at 2 mg/ml (see [Section 3.4.2](#)) are kept at 4°C. Once the vial is opened, its content is distributed in 60-100 µl aliquots in 0.5 ml Eppendorf tubes, stored at -20°C.

Upon application, an equal volume of water is added to one aliquot (in order to prepare 1 mg/ml BSA).

Each aliquot is enough for the positive controls (in the microplates containing the mosquito samples) and for the standard curve, according to alternative 1 (60 µl aliquot) or alternative 2 (100 µl aliquot). *The remaining 1 mg/ml BSA solution is discarded.*

In the standard curve, the point corresponding to 15 µg BSA (see [Sections 5.7.1 and 5.7.2](#)) is obtained directly with the solution at 2 mg/ml. For this purpose, we use a separate aliquot that is re-frozen 2-3 times.

2) **protein reactive (Bio-Rad reactive)** – Bio-Rad cat 500-0006

This solution must be prepared immediately prior to usage.

Working solution

- 8 ml reactive
- 32 ml water

²⁶ We keep several 0.0615 g GSH and 0.0042g CDNB aliquots in Eppendorf tubes, at 4°C.

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(usamos fosfato anidro: Na_2HPO_4)

142 g de fosfato de sódio bibásico

7. GENERAL BUFFERS AND STOCK SOLUTIONS

All the solutions described in **Sections 7.1 to 7.4** are stored at 4°C.

Important notes:

- 1) The buffer stock solution should be at most 20 times more concentrated than the working solution.
- 2) To obtain the required pH, two solutions of identical molarity (mono and bibasic) have to be mixed, in the indicated quantities. The resulting solution (buffer) will have, automatically, the same molarity than the original phosphate solution.
- 3) Mix both solutions in the pHmeter. The amounts indicated below are just an indication.
 - Add more bibasic phosphate if you have to increase pH.
 - Add more monobasic phosphate if you have to decrease pH.
- 4) *Take care with contaminations.* Store your solutions in transparent flasks. Examine always before using!!

7.1. Sodium phosphate buffers

Sodium phosphate buffers needed:

molarity	pH	enzyme, test	solution
100 mM	7.0	ACE	DTNB/Na phosph
50 mM	7.4	PNPA-EST	PNPA
100 mM	7.8	ACE	Triton/Na phosph
20 mM	7.2	EST (alpha, beta)	naphthyl acet/Na phosph

Preparing stocks of sodium phosphate buffers:

7.1.1. sodium phosphate stocks (mono and bibasic) at 1 M

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water to 1 liter

Sterilize by filtration through a 0.22 µm filter

7.1.2. sodium phosphate stock buffers at 1 M

a) The table below discriminates the amounts recommended to prepare 1 liter of sodium phosphate buffers at defined pHs. It is important to assure that the concentration of both phosphate solutions is the same:

pH	Monobasic (NaH ₂ PO ₄)	bibasic (Na ₂ HPO ₄)
7.0	390 ml	610 ml
7.2	280 ml	720 ml
7.4	190 ml	810 ml
7.8	85 ml	915 ml

b) The quantities of 1 M sodium phosphate we use in laboratory routine are as follows:

pH	monobasic (NaH ₂ PO ₄)	bibasic (Na ₂ HPO ₄)
7.0	19.5 ml	30.5 ml
7.2	70 ml	180 ml
7.4	9.5 ml	40.5 ml
7.8	42.5 ml	457.5 ml

c) Use 1 M buffer solutions to prepare buffers that are at least 50 mM.

7.1.3. sodium phosphate buffers at 100 and 50 mM

test	mM	pH	ml buffer 1 M (*)	ml water	ml total
ACE	100	7.0	50	450	500
PNPA-EST	50	7.4	50	950	1000

(*) 1 M buffer solution (Section 7.1.2) in the indicated pH

7.1.4. sodium phosphate buffer at 20 mM

This buffer is used for the alpha and beta-Esterases tests.

To prepare 20 mM sodium phosphate buffer pH 7.2, use stocks of 0.2 M sodium phosphate (mono and bibasic):

a) Dilute the 1 M sodium phosphate stocks 5 times:

- monobasic phosphate: 20 ml phosphate + 80 ml water (100 ml final)

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228,2 g de fosfato de potássio bibásico

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Esterilize com filtro 0,22 μm

- bibasic phosphate: 40 ml phosphate + 160 ml water (200 ml final)
(these will be 0.2 M phosphate solutions)
- b) Prepare 250 ml of 0.2 M buffer, using the same proportions indicated in the table of **Section 7.1.2, b**:
 - 70 ml monobasic phosphate
 - 180 ml bibasic phosphate
- c) Correct the pH, if necessary, using the 0.2 M mono or bibasic sodium phosphate solutions (see **item a**).
- d) Discard the remaining 0.2 M mono and bibasic phosphate solutions (described in **item a**).
- e) Dilute the 0.2 M buffer 10 times to prepare the 20 mM solution:
 - 10 ml 0.2 M phosphate buffer pH 7.2 (see **item c**)
 - 90 ml water
 - This buffer can be diluted in amounts sufficient for usage during 2-3 days.
 - Depending on the work volume, 200 ml can be prepared, mixing 20 ml 0.2 M phosphate buffer and 180 ml water.

7.2. Potassium phosphate buffers

Potassium phosphate buffers needed:

molarity	pH	enzyme, test	solution
90 mM	7.2	MFO	K phosph
100 mM	6.5	GST	GSH

Preparing stocks of potassium phosphate buffers:

7.2.1. potassium phosphate stocks (mono and bibasic) at 1M

- potassium phosphate monobasic 1M – 1 liter (MW = 136.1)

(we use anhydrous phosphate: KH_2PO_4)

136.1 g potassium phosphate monobasic

water to 1 liter

Sterilize by filtration through a 0.22 μm filter

- potassium phosphate bibasic 1 M – 1 liter (MW = 228.2)

(we use trihydrated phosphate: $\text{K}_2\text{HPO}_4 \cdot 3\text{H}_2\text{O}$)

228.2 g de potassium phosphate bibasic

water to 1 liter

Sterilize by filtration through a 0.22 μm filter

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7.2.2. potassium phosphate stock buffers at 1 M

a) The table below discriminates the quantities recommended to prepare 1 liter of potassium phosphate buffers at defined pHs. It is important to assure that the concentration of both phosphate solutions is the same:

pH	monobasic (KH_2PO_4)	bibasic (K_2HPO_4)
6.5	640 ml	360 ml
7.2	283 ml	717 ml

7.2.3. potassium phosphate buffer at 90 mM pH 7.2

a) preparation of buffer at 1 M

- 71.7 ml 1 M bibasic phosphate
- 28.3 ml 1 M monobasic phosphate
- mix on the magnetic stirrer, accompanying with the pHmeter
- correct the pH, if necessary, adding bibasic phosphate to increase the pH or monobasic phosphate, to decrease the pH
- this corresponds to 100 ml of a 1 M potassium phosphate buffer pH 7.2 solution

b) preparation of buffer at 90 mM

- 45 ml 1 M potassium phosphate buffer pH 7.2 (see [item a](#))
- 455 ml water
- check pH. Correct if necessary, using the 90 mM mono or bibasic phosphate solutions, to keep the desired molarity²⁷

7.2.4. potassium phosphate buffer at 100 mM pH 6.5

a) preparation of buffer at 1 M

- 18 ml 1 M bibasic phosphate
- 32 ml 1 M monobasic phosphate
- mix on the magnetic stirrer, accompanying with the pHmeter
- correct the pH, if necessary, adding bibasic phosphate to increase the pH or monobasic phosphate, to decrease the pH
- you will have 50 ml of a 1 M potassium phosphate buffer pH 6.5

b) preparation of buffer at 100 mM

- 50 ml 1 M potassium phosphate buffer pH 6.5
- 450 ml water
- check pH. Correct if necessary, using the 100 mM mono or bibasic phosphate solutions, to keep the desired molarity²⁸

²⁷ Prepare mono or bibasic potassium phosphate at 90 mM mixing 9 ml of 1 M stocks ([Section 7.2.1](#)) and 91 ml water.

²⁸ To prepare 100 mM potassium phosphate mono or bibasic mix 10 ml of the 1 M stocks ([Section 7.2.1](#)) to 90 ml water.

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7.3.1. sodium acetate buffer, stock at 250 mM pH 5.0

a) preparation of buffer at 3 M:

- weigh 246.09 g sodium acetate, anhydrous (MW 82.03)
- dissolve in 800 ml water
- adjust pH to 5.0 with glacial acetic acid
- add water up to 1 liter
- Sterilize by filtration through a 0.22 μm filter

b) preparation of buffer at 250 mM

- 41.6 ml 3 M sodium acetate buffer
- add water up to 450 ml
- check pH again. Correct with glacial acetic acid
- complete volume to 500 ml with water

7.4. Propoxur (Acetylcholinesterase assay)

0.1 M propoxur in acetone

Stock, 1 ml:

- 0.02092 g propoxur (MW: 209.24)
- 1 ml acetone
- prepare 100 μl aliquots and store at 4°C
(6 μl are used in each test)

7.5. SDS (alpha and beta-Esterases test)

5% SDS

In the laboratory this solution is obtained by dilution of the 15% stock (50 ml 10% SDS and 50 ml water)

Stock at 10%, 500 ml:

- 50 g SDS (use mask when weighing!)
- water up to 500 ml
- dissolve with the aid of a magnetic stirrer, *slowly* (to avoid bubbles)
- Sterilize by filtration through a 0.22 μm filter (*also slowly*)
- store at room temperature
(SDS precipitates if stored at 4°C)

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(*) K: leitura cinética; EP: leitura *end point*

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8. COMPARING WHO – CDC – LAFICAVE PROTOCOLS

WHO METHOD (200UL HOMOGENATE)

ENZYME	sample vol (µl)	n° replicate	λ (Abs)	t lecture (minutes)	K X EP (*)	order (**)
EST alpha	20	1	570	15	EP	not indicated
EST beta	20	1	570	15	EP	
PNPA-EST	10	2	405	2	K	
GST	10	2	340	5 K/20 EP	EP/K	
MFO	2	2	650	120	EP	
ACE	25 (x 2)	1	405	5 K/60 EP	K/EP	
PTN	10	2	620	5	EP	

CDC METHOD (2000ul homogenate)

ENZYME	sample vol (µl)	n° replicate	λ (Abs)	t lecture (minutes)	K X EP (*)	order (**)
EST alpha	100	3	540	20	EP	1
EST beta	100	3	540	20	EP	1
PNPA-EST	test not done					
GST	100	3	340	0 – 5	EP	4
MFO	100	3	620	10	EP	3
ACE	100	3 x 2	414	20	EP	2
PTN	20	3	620	0	EP	5


LAFICAVE METHOD (homogenate 300ul)

ENZYME	sample vol (µl)	n° replicate	λ (Abs)	t lecture (minutes)	K X EP (*)	order (**)
EST alpha	10	2	570	15	EP	4
EST beta	10	2	570	15	EP	3
PNPA-EST	10	2	405	2	K	5
GST	15	2	340	20	K	7
MFO	20	2	650	90	EP	1
ACE	25 (x 2)	2	405	60	EP	2
PTN	10	2	620	5	EP	6

(*) K: reading of absorbance as a rate, kinetic lecture; EP: "end point" reading of absorbance

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






ROCK

PTN STANDARD CURVE

Resulting coloration of each assay as indicated above the microplates. PNPA-EST and GST reactions are presented by the same microplate since no color (GST) or a weak yellow color (PNPA-EST) is developed. In the panel "Standard Protein Curve", both alternatives indicated in the text are presented. B: blank; C: positive control; ROCK: susceptible strain (Rockefeller)

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9. ANALYSIS OF RESULTS

At the end of each assay the absorbance value of the replicates of each mosquito is obtained. In order to express these values as enzymatic activity the data must be processed and corrected according to:

- 1) the volume of mosquito homogenization
- 2) the amount of total proteins in each mosquito
- 3) the activity unit considered for each enzyme

In the following items the units used to represent results with each enzyme are indicated. The steps needed to convert absorbance values in enzymatic activity for each mosquito are also briefly discriminated.

9.1. Acetylcholinesterase (ACE)

Acetylcholinesterase results are expressed as the percentage of remaining activity after addition of the inhibitor. The difference between activity in the presence (plate AChI, [Section 5.1](#)) and absence (plate AChE, [Section 5.1](#)) of the inhibitor, propoxur, is calculated:

$$\% \text{ remaining activity} = \frac{\text{AChI} \times 100}{\text{AChE}}$$

Alternatively the result can be expressed as the Acetylcholinesterase inhibition rate:

$$\% \text{ inhibition} = \frac{(\text{AChE} - \text{AChI}) \times 100}{\text{AChE}}$$

In the case of ACE it is unnecessary to correct values according to the total quantity of proteins or to the homogenization volume of each mosquito.

9.2. Total proteins (PTN)

Dosage of total proteins of each mosquito is a necessary step to correct activity values obtained for all enzymes related to metabolic resistance.

To calculate total protein quantity in each mosquito, the absorbance values obtained in the 10 µl aliquots should:

- 1) be corrected by the total homogenate volume (300 µl) –absorbance multiplied by a factor of 30;
- 2) be transformed in µg protein. *In order to do this, it is necessary to divide by the conversion factor.* This factor is obtained with the standard curve ([Sections 5.7.1](#)

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
and 5.7.2), that is used to calculate a linear equation ($y = ax + b$)²⁹. The conversion factor (a) is the slope of this line, and it is obtained by the formula:

$$a = \frac{\Delta y - b}{\Delta x}$$

Where:


Δy = difference between the absorbance values of two points

Δx = difference between the protein values (BSA) of the these points

b = absorbance value in the absence of BSA (blank)

9.3. Mixed Function Oxidases (MFO)

This assay measures the heme content of the insect. In non blood-fed mosquitoes, heme is mainly associated with cytochrome P450 (Hemingway 1998). Cytochrome P450 and its prosthetic heme group are involved in monooxygenation reactions, which catalyze the insertion of an oxygen atom of O₂ in the substrate; the other oxygen atom is reduced to water. The general equation of this reaction is:

The cytochrome P450 enzymes are also known as Monooxygenases cytochrome P450, Mixed Function Oxidases (MFOs), poly-substrate Monooxygenases (PSMOs) and “heme thiolate proteins” or simply P450 (Feyereisen 1999). The assay described in Section 5.2 does not measure the monooxygenation activity of P450 over a substrate, directly. Instead, it quantifies the total heme content of the mosquito, and this enables an indirect estimate of the enzyme activity. Hence, the higher the heme content, higher should be the P450 activity. It is important to note that non blood-fed mosquitoes must be used in order to reduce the number of heme molecules (like those present in the hemoglobin protein) that can interfere with results.

A standard curve is constructed (Section 5.2.1) with cytochrome, containing heme. This standard curve enables conversion of the absorbance values of mosquito homogenates to cytochrome content. Results of this assay are represented as cytochrome μg or nmoles present per total protein of each mosquito, in milligrams ($\mu\text{g cit/mg ptn}$ or nmoles cit/mg ptn).

To calculate the cytochrome concentration/protein mg in each mosquito, the absorbance values obtained in the 20 μl aliquots should:

- 1) be corrected by the total homogenate volume (300 μl) –absorbance multiplied by a factor of 15;
- 2) be transformed in cytochrome μg (or nmoles). In order to do this it is necessary to divide by the conversion factor. This factor is obtained with the standard curve (Section 5.2.1), and calculations are done exactly as described in the previous section;

²⁹ In this equation, y = absorbance, x = BSA concentration in the standard curve, in micrograms, b = absorbance in the absence of protein (“blank”).

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- 3) be corrected by the total protein content of each mosquito. It is then necessary to divide the value obtained above by the total protein content of the corresponding mosquito (that is described above, [Section 9.2](#)).

9.4. alpha and beta-Esterases (α and β -EST)

Results obtained with these enzymes are represented as nmol/mg ptn/min. The absorbance values obtained with 10 μ l mosquito aliquots should:

- 1) be corrected by the total homogenate volume (300 μ l) absorbance multiplied by a factor of 30;
- 2) be transformed in substrate nmols. In order to do this it is necessary to divide by the conversion factors. These factors are obtained with the alpha or beta-naphthol standard curves ([Sections 5.3.1 and 5.4.1](#)). Calculations are done exactly as described in [Section 9.2](#);
- 3) be corrected by total proteins. Each value obtained above is divided by the total protein content of the corresponding mosquito, detailed in [Section 9.2](#). The result is the amount of substrate (in nmol) consumed by protein mg of each mosquito in the total reaction time (15 minutes);
- 4) to calculate the enzyme activity per minute, it is sufficient to divide the value obtained above by 15.

9.5. "PNPA" Esterase (PNPA-EST)

In opposition to alpha- and beta-Esterases, for this assay ([Section 5.6](#)) a standard substrate is not available. For this reason it is neither possible to construct a standard curve nor to convert absorbance values in consumed substrate amount. In this assay the final result is expressed as Δ Abs/mg ptn/min, meaning that the absorbance variation in one minute is calculated, being corrected by the total protein content in each mosquito. In practice the unit that is considered is the rate of reaction (expressed as Δ Abs/min).

The procedure we adopt with the absorbance values obtained in the 10 μ l aliquots is as follows:

- 1) for each mosquito a graph "absorbance X minute" is constructed. With these graphs, the linear equation is calculated, and a correlation index (R^2) is obtained. This index, which varies from 0 to 1, evaluates the correlation between the observed and expected values;
- 2) all mosquitoes that have R^2 values lower than 0.9 are discarded;
- 3) the absorbance difference from 90 to 30 seconds is calculated. This corresponds to the one minute sample variation of reaction (60");
- 4) to obtain the Δ Abs/min for each mosquito, the value obtained above is multiplied by 30 (the dilution factor);

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5) the correction by total proteins is then performed. Each value obtained above is divided by the total protein content of the corresponding mosquito, in milligrams.

9.6. Glutathione-S-transferase (GST)

Similarly to “PNPA” Esterase, for this assay (Section 5.6) a standard substrate is not available. For this reason it is not possible to construct a standard curve. However, in this case it is possible to calculate the amount of substrate used in the reaction. The parameters utilized by Hemingway (1998) have been considered and were incorporated in the calculations. These parameters include:

- the extinction coefficient value of the reaction product, (3-(2-chloro-4-nitrophenyl)-glutathione);
- the spectrophotometer path length, that corresponds to the height of the solution in the microplate well (in our case, 0.6 cm).

These parameters are used to transform absorbance values (Abs) in mmoles/mosquito/minute, according to the formula:

$$\text{mmol/mosq/min} = \frac{(\text{Abs}_{20} - \text{Abs}_{10}) \times 20}{4.39 \times 0.6 \times 10}$$

Where:

Abs_{20} = Abs value obtained after 20 minutes of reaction

Abs_{10} = Abs value obtained after 10 minutes of reaction³⁰

20 in the dividend = dilution factor (15 μl sample, from 300 μl homogenate)

10 in the divisor = minutes of reaction

4.39 = extinction coefficient of the reaction product

0.6 = path length (= height of the solution in the microplate well, in cm)

Simplifying:

$$\text{mmol/mosq/min} = \frac{(\text{Abs}_{20} - \text{Abs}_{10}) \times 2}{4.39 \times 0.6}$$

In order to correct by total proteins, the value above is divided by the total protein content of the corresponding mosquito, obtained as described in Section 9.2. After correction by protein mg, the final result of this assay is expressed as mmols/mg ptn/min, corresponding to the amount, in mmols, of the reaction product generated per minute, per protein mg.

³⁰ In our experience, a high oscillation is observed in the values measured up to 10 minutes of reaction. For this reason we consider only the interval between 10 and 20 minutes.

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10. LAFICAVE METHODOLOGY FOR THE PROCESSING AND ANALYSIS OF RESULTS

The assays herein described generate a great volume of data, since:

- the activity of six enzymes is simultaneously evaluated in each mosquito;
- 40 mosquitoes from the population under test are processed in each assay;
- 80-120 *Aedes aegypti* mosquitoes of each population are evaluated.

To expedite processing of the data generated by the spectrophotometer, four complementary Excel files were elaborated. Together these files enable the transformation of the original data (expressed in absorbance values) in enzymatic activity and also the correction according to the total protein content of each mosquito. In addition, these files facilitate the calculations discriminated in the previous sections. Their final result is the automatic construction of histograms – one for each enzyme of the populations evaluated. In these histograms the *x*-coordinate refers to the activity of the enzyme (shown in classes), and the *y*-coordinate exhibits the number of specimens in each class. In this way, at the end of each analysis it is possible to evaluate the enzymatic profile of a given population, which may then be compared to the profile of the susceptible reference strain (Figure 1).

In the following sections, the general procedures for usage of each Excel file are described³¹. Detailed instructions for filling out **Files 1 to 4** are included in the first sheet (*instructions*) of each file.


Figure 1: typical histogram generated by **File 2**. alpha-Esterase activity profile of the susceptible strain (upper graph) and an altered field population (lower graph).

³¹ These files are available in the CD that accompanies the present handbook.

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10.1. File 1 – data transfer from Soft Max to Excel

The spectrophotometer we use in the laboratory is associated with the Soft Max program that shows the absorbance value for each microplate position. However, as mentioned above, the files for data processing were elaborated in the Excel program.

When data are transferred directly from Soft Max to Excel, the absorbance values are accompanied by other parameters, such as temperature and reading time, which are not necessary for the processing steps. These “excess” data should be eliminated, and only absorbance data should be retained.

The principal purpose of **File 1** is to prepare original data in such a way only the absorbance values are maintained. This procedure facilitates transference of data to the other files. By doing this it is also possible to work in other computers, where Soft Max is not installed.

File 1 contains six sheets. The first one presents filling out instructions; each of the remaining sheets was designed to receive data from all the enzymes of a complete test. These sheets are referred to as “*date1*” to “*date5*” and should be renamed with the date of the assay that will be inserted.

Each sheet has cells assigned to receive the original absorbance values from all the enzymes, from the quantification of total proteins and from the protein standard curve (BSA). These cells, that are blue-labeled and simulate the 96-well microplates, follow the numbering and distribution of the mosquitoes previously defined: mosquitoes 1 to 40 correspond to the population under test, and mosquitoes 41 to 45 belong to the susceptible reference strain, as described in **Section 4.3.2**.

How to use File 1:

- 1) according to the sheet “*instructions*” of **File 1**, save as the file, with the expression “*original data*”, followed by the name and localization of the analyzed population and the year of collection of specimens;
- 2) rename each sheet “*date*” with the date of the test: month, day, year;
- 3) copy the original data (provided by the spectrophotometer) from the Soft Max program and paste in blank cells on the sheet. Then select only cells containing the absorbance values. These values are transferred to the blue-labeled cells, designated to receive the absorbance values for each enzyme.

Henceforth, data inserted in **File 1** sheets will be used to fill out the following files (**Files 2 to 4**). These files contain calculations to convert the absorbance values in enzymatic activity (**Files 2 and 3**) and to pool all the tests from the same population (**File 4**).

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SoftMax Pro - [PNPA_25.09.2006.pda]

File Edit View Experiment Control Assays Data Window Help

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


10.2. File 2 – conversion of absorbance values (except PNPA-EST)

In this file original data are converted to enzymatic activity and corrected according to the total protein content of each mosquito. This file includes all the calculations described in [Section 9](#) for each tested enzyme with the exception of PNPA-EST, calculated in a separate file ([File 3](#)).

[File 2](#) has 14 sheets, the first of which, “instructions”, contains the filling out instructions. The 10 following sheets are designated for the data processing of each specific assay. These sheets correspond to the enzyme evaluated, in some cases more than one sheet for the same assay:

- protein standard curve – two sheets, if the standard curve, with BSA, has been established according to alternative 1 or 2 ([Sections 5.7.1 or 5.7.2](#));
- GST – two sheets: the first, “GST kinetics”, for transference of all the data from this assay and the other, “GST T20-T10”, for the timepoints that are used in the calculations (10 and 20 minutes of reaction). This second sheet is completed automatically, when data are transferred from [File 1](#) to [File 2](#), sheet “GST kinetics”;
- MFO – two sheets: the first for activity calculation when expressed in μg (“MFO μgC ”) and the other in nmols (“MFO nmolC”).

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Three additional sheets follow:

- a) “*consolidate*”, a sheet that automatically exhibits (with exception of PNPA-EST, see [Section 10.3](#)) a table with the final results for all mosquitoes analyzed for each enzyme that has been quantified;
- b) “*graphs*”, a sheet that presents tables (one for enzyme) with the classification of the evaluated mosquitoes with respect to range of activity. This sheet also presents histograms, one for each enzyme, with the profiles of the tested mosquitoes. These histograms display the profiles of the population evaluated (40 mosquitoes per test, maximum) and mosquito controls (5 mosquitoes, maximum);
- c) “*conventions*”, this sheet contains information regarding the mosquito homogenate volumes used in the assays with each enzyme and provides instructions for the construction of the histograms with the population profiles of the enzymatic activities. This sheet is not routinely used and was included only for consultation.

Similar to [File 1](#), each sheet of this file has blue-labeled areas, simulating the 96-well microplate, where [File 1](#) values are “pasted”. The positions corresponding to reference values (“blanks”, positive controls and reference strain mosquitoes, respectively in white, yellow and pink) are also discriminated. In some [File 2](#) sheets there are lilac cells, corresponding to the places where the slopes of some standard curves must be inserted, as detailed below (“*Observations regarding File 2*”, *items 1 and 3*).

How to use File 2:

- 1) according to the “*instructions*” sheet of [File 2](#), save as the file, with the name and location of the analyzed population, together with the date of the test;
- 2) transfer data from [File 1](#) to the sheets corresponding to each enzyme of [File 2](#). This is done for all the enzymes with the exception of PNPA-EST (see [Section 10.3](#));
- 3) discard samples if replicates present standard deviation beyond the acceptable. Follow indications on the “*instructions*” sheet;
- 4) verify the profiles obtained on the sheet “*graphs*”. Data for all the enzymes are also displayed in one table on the sheet “*consolidate*”. Data in this table can be transferred to several other programs if one wants to further analyze or construct other kinds of graphs.
- 5) The criteria utilized in our laboratory for the interpretation of results are detailed in [Section 11](#).

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
Observations regarding File 2:

- 1) one protein standard curve, with BSA, is established for each new biochemical assay. This is the reason for a “*ptn curve*” sheet in **File 2**. Data obtained with this curve are used in the dosage of total proteins, and the slope (a) value of the standard curve is utilized to calculate the amount of total proteins of each mosquito of a given test. From this point on all the enzymatic activity values will be automatically corrected according to the total protein content of each mosquito in all the enzyme sheets related to metabolic resistance;
- 2) only in the case of Acetylcholinesterase, the enzyme related to target-site resistance, is correction according to total protein content not necessary: analysis of these results is based on the difference in activity with or without propoxur (**Section 9.1**);
- 3) the slope values obtained with the Esterases and MFO standard curves (**Sections 5.3.1, 5.4.1 and 5.2.1**) are also used to convert absorbance in activity values in the corresponding sheets of **File 2**. Note that these standard curves are not included in this file, as they are constructed separately. To calculate the slopes of these curves, the same procedure for the protein standard curve (described in **Section 9.2**) is adopted. The construction frequency of these curves depends on individual laboratory necessity (see *footnote 14*, at **Section 5**);
- 4) in the lower part of the enzyme sheets, the final result related to each mosquito, already corrected, will be displayed. These values will be repeated in the “*consolidate*” sheet, that pools results of all **File 2** sheets;
- 5) the sheet “*graphs*” exhibits the profile of each enzyme, for control and tested mosquitoes, thereby enabling comparison (see **Section 11** for details);
- 6) it is important to stress that organization of all the results in one table allows the analysis of data in different ways, as well as affording representations different from the histogram displayed in **Figure 1**. For instance, the consolidation of these data enables same mosquito evaluation of activity related to each enzyme and establishment of correlations among the enzymes that may be relevant. In addition, other graphic representations can be adopted (see **Section 11** for examples);

Microsoft Excel - TBL 2 - Conversion absorption values

The screenshot shows a Microsoft Excel spreadsheet titled "TBL 2 - Conversion absorption values". The spreadsheet contains a large table with multiple columns and rows. The columns are labeled "Absorption" and "Conversion", with sub-columns for "100%", "50%", and "25%". The rows are numbered 1 through 30. A red circle highlights a cell in the bottom right corner of the table, containing the text "comodora".

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10.3. File 3 – conversion of PNPA-EST absorbance values

We consider that, in the case of PNPA-EST, it is important to include only mosquitoes whose activity kinetics does not oscillate excessively. To accomplish this, it is necessary to verify, for every mosquito, the correlation between the obtained results and the expected values. This is achieved by the construction of graphics exhibiting the enzymatic activity of each mosquito throughout time together with the calculation of the respective correlation indexes (R^2). Since this evaluation generates a huge volume of data, the conversion of PNPA-EST absorbance values is reserved for a separate file.

File 3 calculates, for all the mosquitoes of a given test, the activity of PNPA-EST, expressed in $\Delta\text{Abs}/\text{mg protein}/\text{min}$.



File 3 contains 3 sheets, entitled:

- a) “*instructions*”;
- b) “*PNPA calculation*”, containing graphics for individual mosquitoes. In this sheet it is possible to evaluate the mosquitoes presenting PNPA-EST kinetics with a high oscillation. It is also possible to discriminate samples with an unacceptable standard deviation;
- c) “*activity calculation*” where correction of activity according to the total protein amount of each mosquito is performed. Graphs with the PNPA-EST activity profile of the pool of mosquitoes evaluated in one test is displayed is shown in File 2, sheet “*graphs*”.

The sheet “*PNPA calculation*” presents blue-labeled cells assigned to receive the original absorbance values of the whole kinetics of this enzyme. Total protein values of each mosquito, calculated in **File 2**, are “pasted” in the sheet “*activity calculation*”.

How to use File 3:

- 1) according to the sheet “*instructions*” of **File 3**, save as the file with the expression “PNPA”, followed by the name and location of the analyzed population together with the date of the test;
- 2) transfer PNPA-EST kinetics results from **File 1** to the sheet “*PNPA calculation*”;
- 3) transfer total protein/mosquito data from **File 2** to the sheet “*activity calculation*”;
- 4) discard samples if replicates present standard deviation beyond acceptable or if the correlation index is low. Follow indications in the sheet “*instructions*”;
- 5) copy the results obtained from the sheet “*activity calculation*” and paste to **File 2** of the corresponding test, sheet “*consolidate*”.

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Instructions for filling out the spreadsheet

- 1) open this file and "save as" in the corresponding population folder
The file name will be: "PNPA [population name] [localization] [date of the test: Month Day Year]"
examples:
A) if it is a municipality:
"PNPA Los Angeles CA 01 25 2003"
B) if it is a neighborhood:
"PNPA Hollywood LA CA 01 24 2003"
In the case of the reference strain insert only "PNPA" followed by the strain abbreviation and the date of the test example (for Rockefeller strain): "PNPA Rock 01 25 2006"
- 2) paste PNPA kinetics data (0 to 120 seconds) from **File 4** to the sheet "PNPA calculation"
- 3) evaluate each mosquito in the graphs formed in the "PNPA calculation" sheet:
discard mosquitoes if the correlation index is lower than 0.9
The correlation index value (R^2) is in the upper left border of each mosquito's graph in the sheet "PNPA calculation"
also discard mosquitoes for which the standard deviation is higher than 30% of the mean in 5 or more points of the kinetics
this calculation appears in the tables "relation mean/stddev", starting at cell "AV35", in the sheet "PNPA calculation", where all the values to be discarded will appear in **bold and red**

To discard the mosquito:
in the sheet "activity calculation" eliminate the value of the corresponding mosquito from the "I" column



- 3) insert, in the sheet "activity calculation", the total protein dosage values from the corresponding test copy, from **File 2** of the corresponding test sheet "consolidate": the total protein values (these values are placed in cells "D9" to "D63")
paste special (as "values"), the total protein values, to this **File 3**, sheet "activity calculation", starting at cell "F5"
if any protein value is not acceptable, remove, from the "I" column of this sheet, the corresponding mosquito result
- 4) copy the PNPA activity values to **File 2**, sheet "consolidate"
copy from **File 3**, "activity calculation" sheet, cells "B5" to "M49" to **File 2**, "consolidate" sheet, cells "B5" to "M43"

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



10.4. File 4 – total data of a given population



File 4 pools the results of all the assays with one population. In principle, one **File 4** should be created for each population. Alternatively, there may be only one **File 4** containing several sheets, one for each population evaluated in a given period of time.


This file consists of, besides the sheet “*instructions*”, two more sheets: “*data*” and “*graphs*”. Results obtained from all the completed tests (**File 2**, after insertion of PNPA-EST) are transferred to the sheet “*data*”. The sheet “*graphs*” generates, automatically, histograms with the profiles of the population for each enzyme evaluated.

How to use File 4:

- 1) according to the sheet “*instructions*” of **File 4**, “*save as*” the file with the expression “total”, followed by the name and location of the analyzed population with the year of collection of specimens;
- 2) alternatively, a “*total*” file can be generated for a group of populations evaluated in a given period. In this case, the file should be entitled “total” followed by the year of specimen collection;
- 3) the sheet “*data*” has space, organized in columns, to insert results obtained with each enzyme of all the tests with the same population. In order to do this it is sufficient to paste in these cells data “*copied*” from the sheets “*consolidate*” of each **File 2** of that population: data concerning the 40 mosquitoes from the population under test (numbers 1 to 40, since mosquitoes 41 to 45 belong to the control strain);
- 4) this procedure is repeated, sequentially, with every test performed, until this sheet contains the results from all the tests with a given population;
- 5) the enzymatic activity results and the number of tested mosquitoes are used in the construction of histograms that appear simultaneously with the insertion of data, on the sheet “*graphs*”. As a final result, the enzymatic profile of the population can be observed;
- 6) enzymatic activity classes defined in **File 4** are adequate for the mosquito populations that are being evaluated in our laboratory. If the profile of a given population is better represented by other activity interval classes, it is sufficient to redefine them. In order to do this, go to **File 2**, sheet “*conventions*”.

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
11. INTERPRETATION OF RESULTS

11.1. Criteria for inclusion of evaluations from individual mosquitoes

- 1) In order to control the quality of manipulation, all the assays are carried out with duplicate samples. Accordingly, we decided as an arbitrary criterion, to reject all samples for which the standard deviation between the replicas was higher than 30% of the mean. In **Files 2** and **3** the mean and standard deviation between replicas are calculated automatically. These files also calculate the standard deviation representation with respect to the mean. In addition, samples that are to be rejected appear in bold and red. Suffice it to eliminate these samples from the cells indicated in **Sections 10.2** (“How to use **File 2**”, item 3) and **10.3** (“How to use **File 3**”, item 4) as well as on the “instructions” sheet of these files.
- 2) Quantification of PNPA-EST evaluates the reaction kinetics of the enzyme. In this case, in addition to the rejection of samples according to the criterion detailed above, we opted to eliminate mosquitoes when the oscillation of their reaction kinetics is very high. We defined, also as an arbitrary criterion, the elimination of samples when the correlation between the obtained and expected results (R^2) is lower than 0.9. Filling out **File 3** with absorbance data generates one graph for

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each evaluated mosquito. In these graphics, the R^2 value must be inspected and samples with R^2 lower than 0.9 are discarded from the cells indicated in **Section 10.3** (“How to use **File 3**”, item 4) and on the “instructions” sheet.

- 3) Prior to test initiation with field populations it is necessary to obtain profiles for all the enzymes herein discriminated of an insecticide-susceptible reference strain. This profile is obtained after at least three complete assays (totalizing approximately 100-120 mosquitoes). The profiles of the reference strain for each enzyme enable comparisons with field populations (our laboratory uses Rockefeller strain).
- 4) The reference strain profiles mentioned above also help to control the quality of the tests performed with field populations. To accomplish this, it is convenient to compare the reference strain profiles, mentioned above, with the results obtained from the 5 mosquitoes used as internal controls in each test. If values are equivalent, the test can be considered valid.

11.2. Criteria for population classification

1) **Graph representation:** one option to the conventional graphic representation, shown in **Figure 2**, is the inclusion, in only one graph, of all populations from a series of evaluations. This procedure facilitates the direct comparison of enzymatic activity among populations. In this graph, the inclusion of the insecticide-susceptible reference strain enables visualization of the activity alteration level of each population.

However we consider that for this kind of representation parameters such as mean and standard deviation are not adequate: in several cases individuals of one population do not follow the normal distribution for the characteristics evaluated in the scope of this handbook. In particular, populations with an incipient enzymatic alteration, presenting only a few affected specimens, are especially relevant to resistance status monitoring. Identification of these populations allows intervention in the field before resistance attains a threshold that represents a risk to vector control.

One of the alternatives to the parameters “mean” and “standard deviation” is the use of the median. In this case, different populations are represented in box-plot graphs, where the line inside the box indicates the median, the box limits represent the quartiles and the vertical lines outside the box exhibit the limits of the distribution (**Figure 2 a**). Another graph representation that may also be adopted appears as a cloud of points (scatter plot). This representation displays directly the activity of each individual, giving a precise description of the population profile (**Figure 2 b**).

Neither of these two representations is available in **Files 1 to 4**, herein described. Our laboratory procedure is to transfer data from **File 4** to program GraphPad Prism version 4.00, in order to construct the box or scatter plot.

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
Figure 2: examples of graph representation of enzymatic activity in *Aedes aegypti* populations. Panel on top (a): box-plot; panel above (b): scatter plot. “PNPA” Esterase activity profiles are shown in the reference strain, Rockefeller (ROCK) and in three field populations (A, B, C).

2) How to classify a profile as normal or altered

Enzymatic profile comparisons among populations through parametric and non-parametric statistical analysis revealed significant differences even when these differences were not functionally relevant. This was, for instance, the case with ACE activity in Brazilian vector populations: statistical analysis demonstrated that several field populations were significantly different from the susceptible reference strain. However, according to WHO criterion (a functional criterion), this enzyme is considered susceptible if the inhibition of its activity in the presence of propoxur is 70% or more (Hemingway, 1998). When this criterion was employed in the evaluated populations, it was verified that ACE activity of all of them should be considered susceptible.

We then decided to work in accordance with the definition of a classification that was representative of these functional differences. Precisely due to this, we first evaluated ACE activity, as this enzyme is associated with a functional classification criterion that is widely accepted (Hemingway, 1998).

A classification was defined that represents adequately the functional differences that were found, regarding not only ACE but also other enzymes. This classification is based on the use of reference strain 99 percentile as the cutoff value. In each population, the percentage of specimens with activity above the 99 percentile of the reference strain is calculated. Some examples are exhibited below:

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Step 1 – calculate, for each enzyme, 99 percentile of the reference strain (this is easily achieved with several spreadsheet softwares, like Excel).

Examples of 99 percentile for a reference strain.

ACE	MFO	GST	α -Est	β -Est	PNPA-Est
(% inhibition activity)	nmoles cit/ mg ptn	mmoles/ mg ptn/min	nmol/ mg ptn/min	nmol/ mg ptn/min	Δ Abs/ mg ptn/min
25,12	33,77	1,79	48,62	81,18	3,58

Step 2 – calculate, for each population, the percentage of specimens with activity higher than the 99 percentile of the reference strain.

The numbers indicate, for hypothetical populations, the percentage of specimens above the 99 percentile of the reference strain.

Population	ACE	MFO	GST	α -EST	β -EST	PNPA-EST
A	0	0	21	14	4	21
B	0	0	45	39	39	44
C	0	0	61	33	2	27
D	2	2	0	23	21	14
E	0	24	0	93	80	0
F	0	0	0	3	1	18
G	3	35	2	77	61	19
H	11	0	15	0	19	81
I	0	0	0	48	16	0

This cutoff value enables the consistent identification of measured activity differences. In addition, it is possible to quantify this difference; since numeric values are assigned to each population (it is relevant to note that these values express a comparison with the reference strain). We then decided, arbitrarily, to classify these values: 0-15 corresponds to “unaltered” activity, equivalent to the reference strain; 15-50 suggest “altered activity”; values above 50 represent “highly altered” activity. It is possible to give colors to these numbers in order to facilitate the result visualization.

População	ACE	Monox	GST	α -EST	β -EST	PNPA-EST
A	0	0	21	14	4	21
B	0	0	45	39	39	44
C	0	0	61	33	2	27
D	2	2	0	23	21	14
E	0	24	0	93	80	0
F	0	0	0	3	1	18
G	3	35	2	77	61	19
H	11	0	15	0	19	81
I	0	0	0	48	16	0

% pop. legenda

0-15 

15-50 

50-100 

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- Feyereisen R 1999. Insect P450 Enzymes. *Annu Rev Entomol.* 44: 507-33.
- Hemingway J 1998. Techniques to detect insecticide resistance mechanisms (field and laboratory manual). WHO/CDC/CPC/MAL/98.6, World Health Organization, Geneva.
- Insecticide Resistance Workshop 1998. Centers for Diseases Control. "Microplate assays", Atlanta, GAS.³²
- Lima-Castelani AR, Ceron CR, Bicudo HE 2004. Variation of genetic expression during development, revealed by esterase patterns in *Aedes aegypti* (Diptera, Culicidae) *Biochem Genet* 42: 69-84.
- Montellano PRO 2004 Cytochrome P450: structure, mechanism and biochemistry. 3rd ed. Springer. ISBN: 0306483246.
- Mourya DT, Hemingway J, Leake CJ 1993. Changes in enzyme titers with age in four geographical strains of *Aedes aegypti* and their association with insecticide resistance. *Med Vet Entomol* 7: 11-16.
- Sousa-Polezzi R, Bicudo H 2005. Genetic variation along time in a Brazilian population of *Aedes aegypti* (Diptera: Culicidae), detected by changes in the esterase patterns. *Genetica* 125: 43-53.

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Brogdon WG, Dickinson CM 1983. A microassay system for measuring Esterase activity and protein concentration in small samples and in high-pressure liquid chromatography eluate fractions. *Analyt Biochem* 131: 499-503.


Brogdon WG 1988. Microassay of Acetylcholinesterase activity in small portions of single mosquito homogenates. *Comp Biochem Physiol* 90C: 145-150.

Brogdon WG, Barber AM 1990. Microplate assay of Glutathione S-transferase activity for resistance detection in single-mosquito triturates. *Comp Biochem Physiol* 96B: 339-342.

Brogdon WG, McAllister JC, Vulule J 1997. Heme Peroxidase activity measured in single mosquitoes identifies individuals expressing an elevated Oxidase for insecticide resistance. *J Amer Mosq Cont Assoc* 13: 233-237.

Population	ACE	MFO	GST	α -EST	β -EST	PNPA-EST
A	0	0	21	14	4	21
B	0	0	45	39	39	44
C	0	0	61	33	2	27
D	2	2	0	23	21	14
E	0	24	0	93	80	0
F	0	0	0	3	1	18
G	3	35	2	77	61	19
H	11	0	15	0	19	81
I	0	0	0	48	16	0

% pop. legend

0-15	
15-50	
50-100	

12. REFERENCES

- Feyereisen R 1999. Insect P450 Enzymes. *Annu Rev Entomol* 44: 507-33.
- Hemingway J 1998. Techniques to detect insecticide resistance mechanisms (field and laboratory manual). WHO/CDC/CPC/MAL/98.6, World Health Organization, Geneva.
- Insecticide Resistance Workshop 1998. Centers for Diseases Control. "Microplate assays", Atlanta, GAS.³²
- Lima-Castelani AR, Ceron CR, Bicudo HE 2004. Variation of genetic expression during development, revealed by esterase patterns in *Aedes aegypti* (Diptera, Culicidae) *Biochem Genet* 42: 69-84.
- Montellano PRO 2004 Cytochrome P450: structure, mechanism and biochemistry. 3rd ed. Springer. ISBN: 0306483246.
- Mourya DT, Hemingway J, Leake CJ 1993. Changes in enzyme titers with age in four geographical strains of *Aedes aegypti* and their association with insecticide resistance. *Med Vet Entomol* 7: 11-16.
- Sousa-Polezzi R, Bicudo H 2005. Genetic variation along time in a Brazilian population of *Aedes aegypti* (Diptera: Culicidae), detected by changes in the esterase patterns. *Genetica* 125: 43-53.

³² In this Workshop slight modifications of the enzymes quantification methodologies, described in the references below, were presented.

Brogdon WG, Dickinson CM 1983. A microassay system for measuring Esterase activity and protein concentration in small samples and in high-pressure liquid chromatography eluate fractions. *Analyt Biochem* 131: 499-503.

Brogdon WG 1988. Microassay of Acetylcholinesterase activity in small portions of single mosquito homogenates. *Comp Biochem Physiol* 90C: 145-150.

Brogdon WG, Barber AM 1990. Microplate assay of Glutathione S-transferase activity for resistance detection in single-mosquito triturates. *Comp Biochem Physiol* 96B: 339-342.

Brogdon WG, McAllister JC, Vulule J 1997. Heme Peroxidase activity measured in single mosquitoes identifies individuals expressing an elevated Oxidase for insecticide resistance. *J Amer Mosq Cont Assoc* 13: 233-237

ISBN 85-334-1291-6

disque saúde
0800.61.1997

www.saude.gov.br/svs

www.saude.gov.br/bvs

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© 2006 Ministério da Saúde.